



CSP Dairy™

3 oz Solid-Top Vial

Single Use Clean Flip-Top Vial for Dairy, Food and Veterinary Sampling

Applications

Fields of use

- Field collection of farm milk and related dairy products
- Quality control sample storage containers and food packaging vials
- Veterinary specimen collection, infection control of mastitis
- Environmental sampling containers

Options

- Can be supplied with label (barcode...)
- Standard vial is colourless - colours upon request

Benefits

Flip-Top Design

- Consumer-friendly concept (easy open/close)
- Patented "in mould closing" process
- Perfectly tight and clean by process
- Translucent in colour for visual inspection
- Pure Polypropylene suitable for food contact
- 100% Recyclable after use

User-Friendly Hinged Vial Concept

- One-hand opening and closing for perfect cleanliness and ease of use
- Allows the use of gloves in cold environment

Regulatory and Quality

- ISO-9001 certified company
- Complies with Regulation (EC) No 1935/2004

Characteristics

Please note that these values are average measurements
Tolerance +/- 0.5 mm (will vary based on product)

Internal volume	77 ml
Filling level line	60 ml
Vial overall height	85,8 mm
Body height	83,5 mm
Inside height	83,7 mm
Top lid diameter	40,1 mm
Top body diameter	36,3 mm
Inside top diameter	35,6 mm
Bottom body diameter	35,1 mm
Vial weight	10,4 grams

Production mould # 723

Packaging

ART. # 03 03 6086 0000 CN

Clear in color
400 Pieces per box
30 Boxes
12.000 Pieces per pallet

