



CSP Dairy™

4 oz High Profile Solid-Top Vial

Single Use Clean Flip-Top Vial for Dairy, Food and Veterinary Sampling

Applications

Fields of use

- Field collection of farm milk and related dairy products
- Quality control sample storage containers and food packaging vials
- Veterinary specimen collection, infection control of mastitis
- Environmental sampling containers

Options

- Can be supplied with label (barcode...)
- Standard vial is colourless - colours upon request

Benefits

Flip-Top Design

- Consumer-friendly concept (easy open/close)
- Patented "in mould closing" process
- Perfectly tight and clean by process
- Translucent in colour for visual inspection
- Pure Polypropylene suitable for food contact
- 100% Recyclable after use

User-Friendly Hinged Vial Concept

- One-hand opening and closing for perfect cleanliness and ease of use
- Allows the use of gloves in cold environment

Regulatory and Quality

- ISO-9001 certified company
- Complies with Regulation (EC) No 1935/2004

Characteristics

Please note that these values are average measurements
Tolerance +/- 0.5 mm (will vary based on product)

Internal volume	131 ml
Filling level line	90 ml
Vial overall height	90,2 mm
Body height	89,2 mm
Inside height	88,1 mm
Top lid diameter	49,5 mm
Top body diameter	45,9 mm
Inside top diameter	44,9 mm
Bottom body diameter	44,9 mm
Vial weight	17,1 grams

Production mould # 800

Packaging

ART. # 03 04 7090 00HP CN

- Clear in color
- 250 Pieces per box
- 30 Boxes
- 7.500 Pieces per pallet



Your product, actively protected

csptechnologies.com +33 (0)3 88 09 89 00 info@csptechnologies.fr

9 rue du Sandholz, 67110 Niederbronn-les-Bains, France

Global Offices Auburn, AL • Amsterdam, NY • Niederbronn-les-Bains, FR • Paris, FR • Osaka, JP • Hong Kong, CN

1. CSP Technologies logo and tagline are service marks of CSP Technologies, Inc., AL.
2. CSP Dairy is trademark of CSP Technologies, Inc., AL.