

TECHNICAL DATA SHEET

| | |
|--------------------|-----------------|
| Part Number | 068865 |
| Description | 3 Fingers clamp |



Benefits :

- ✓ Made from chrome / aluminium
- ✓ Ring for separating funnels and funnels
- ✓ Clamp with removable nut : central adjustment which allows the inclination to be adjusted

Certificate :

- ✓ Non sterile

Technical specifications :

- ✓ Distance : 0-95 mm
- ✓ Rod Ø 8 mm
- ✓ Length 140 mm

Sales unit :

1x 3 Fingers clamp
