



SciLabware

QF/7.17 Issue A 10 Nov 2017

SciLabware Limited  
Unit 4, Riverside 2, Campbell Road  
Stoke-on-Trent, ST4 4RJ, United Kingdom

Tel +44 (0) 1782 444406

Fax +44 (0) 1782 940436

[www.scilabware.com](http://www.scilabware.com)  
[www.wheaton-uk.com](http://www.wheaton-uk.com)

Registered office (as above)  
Registered in England and Wales No. 7723574  
BS EN ISO 9001: 2008 registered company

## PRODUCT SPECIFICATION

<b>Product Code:</b>	1000/28D
<b>Product Description:</b>	2000ml Low Form Beaker
<b>Brand:</b>	PYREX
<b>Type of Glass:</b>	Borosilicate glass Type 3.3, which conforms to ISO3585 European and US Pharmacopoeia Type 1 glass
<b>Product Specification:</b>	White graduations and marking spot  Height 185mm O.D. 130mm
<b>Product Standard(s):</b>	The product is manufactured to comply with ISO3819 and DIN12331
<b>Pack Quantity:</b>	10 pieces per pack

SciLabware Limited hereby confirms that the product specified above is inspected according to all relevant specifications and procedures laid out in our Quality System, in accordance with our registration to BS EN ISO9001 and BS EN ISO13485.

Signed on behalf of SciLabware Limited by:

Signed:   
\_\_\_\_\_  
D Turner  
Quality Technician

Date: 13th Feb 2018

Pyrex® is a registered trademark. SciLabware registered user.

Azlon®

MBL®

PYREX®



SVL®

WHEATON®