

L T C R A N G E



LIQUID

TEMPERATURE CONTROL

Froilabo

SINCE 1918
LYON, FRANCE

Transparent Water Baths

Quality that you can see clearly

- Polycarbonate tank fits for applications required to observe sample changes.
- Compact size that save bench space
- Temperature range from ambient +7°C to 60°C
- Scaled temperature analogue control with commonly use +37°C and +56°C mark.
- Low water level power cut-off
- Optional covers are available for reducing heat loss and preventing evaporation.



Accessories

- Roof style plexiglas cover with thermometer holder for BMP30
- Roof style plexiglas cover with thermometer holder for BMP60
- Holder for 30 hemolysis tubes dia. 14 mm (BMP60 can fit 2 pieces)
- Holder for 30 test tubes dia. 20 mm (BMP60 can fit 2 pieces)

BM Series Specifications

Water Bath Type	Model	Volume (L)	Heater Power (W)	Power Supply (Voltage / Hz)	Exterior Dimension (mm) L x H x D	Interior Dimension (mm) L x H x D
Stainless Steel Analogue	BMUT15	15	1000	220 / 50	390 x 265 x 430	300 x 160 x 280
	BMUT25	25	1500	220 / 50	565 x 265 x 430	490 x 160 x 280
	BMUT45	45	2000	220 / 50	680 x 265 x 630	590 x 160 x 480
Stainless Steel Digital	BMUTE15	15	1000	220 / 50	390 x 265 x 430	300 x 160 x 280
	BMUTE25	25	1500	220 / 50	565 x 265 x 430	490 x 160 x 280
	BMUTE45	45	2000	220 / 50	680 x 265 x 630	590 x 160 x 480
Transparent Analogue	BMP30	3	250	220 / 50	360 x 120 x 130	280 x 110 x 80
	BMP60	6	250	220 / 50	370 x 150 x 210	290 x 110 x 190