

## TECHNICAL DATA SHEET

<b>Part Number</b>	069801
<b>Description</b>	Non-pyrogenic sterile water Otec 500 ml



### Benefits :

- ✓ Aguettant OTEC water is an Active Pharmaceutical Ingredient (API) (Used as a pharmaceutical primary material excluding injectable, ophtalmic or inhaled products)
- ✓ For rinsing and irrigation of materials and tools, white rooms cleaning, machine washing and acting as an active substance for the manufacture of médicaments
- ✓ OTEC range cannot be used neither in the hospital sector, nor for rinsing and irrigation of wounds

### Certificate :

- ✓ Sterile
- ✓ Complies with directives NF EN ISO 3696 type 3 and European pharmacopeia's monograph «purified water in container »

### Technical specifications :

- ✓ Volume : 500 ml

### Sales unit :

1 Bottle