

TECHNICAL DATA SHEET

Référence	069805A
Désignation	Non pyrogenic sterile water



Description :

- ✓ Purified water, sterile and pyrogen-free solution
- ✓ Conductivity at 25°C $\leq 5,0 \mu\text{S}/\text{cm}$
- ✓ Aguettant OTEC water is an Active Pharmaceutical Ingredient (API) (Used as a pharmaceutical primary material excluding injectable, ophtalmic or inhaled products)
- ✓ Complies with directives NF EN ISO 3696 type 3 and European pharmacopeia's monograph «purified water in container»
- ✓ For rinsing and irrigation of materials and tools, white rooms cleaning, machine washing, and acting as an active substance for the manufacture of médicaments
- ✓ OTEC range cannot be used neither in the hospital sector, nor for rinsing and irrigation of wounds.
- ✓ Sterile
- ✓ Volume 5 L
- ✓ Packaged weight 5611 g
- ✓ Packaged dimensions 27 x 21,5 x 20 cm

Sales unit :

1 Bottle