


## Disposable robotic tips 300 µl conductive



### Product description

Conductive tips, 300 µl, tray, 96 pieces in a tray

<b>Part Number</b>	49000-0000
<b>GTIN</b>	N.A.
<b>GDMN</b>	Pipet tip
<b>Tip volume</b>	300 µl
<b>Color Tip</b>	black (conductive)
<b>Color Tray</b>	beige
<b>Sterile</b>	no
<b>Re-usable</b>	no
<b>Patient contact</b>	no
<b>Contact with blood and blood plasma</b>	yes

### Classification

<b>Europe CE labeling</b>	CE-IVD	<i>In accordance with Directive 98/79/EC on in vitro diagnostic medical devices as accessory.</i> <i>Article which, whilst not being an in vitro diagnostic medical device, is intended specifically by Ritter GmbH to be used together with a device to enable that device to be used in accordance with its intended purpose.</i>
---------------------------	--------	--

### Purity



Produced under cleanroom conditions  
(ISO 14644-1 class 8)

## Materials

Tip	PP – Polypropylene
Tray	PC – Polycarbonate


## Tip Dimensions

Tip Length	55.0 mm
Concentricity	≤ 1.5



## Tray Dimensions

Tray length	126.9 mm
Tray width	88.4 mm



## Packaging

Tray	96 tips
Box	180 trays (17,280 tips)
Pallet	12 boxes (207,360 tips)

Specifications are subject to change without prior notice.

In case of any questions, please contact our sales representative via phone, e-mail or our website.

<https://www.ritter-medical.de/en/contact>