


Disposable robotic tips 1100 µl conductive



Product description

Conductive tips, 1100 µl, tray, 96 pieces in a tray

| | |
|--|--------------------|
| Part Number | 49001-0000 |
| GTIN | N.A. |
| GDMN | Pipet tip |
| Tip volume | 1100 µl |
| Color Tip | black (conductive) |
| Color Tray | gray |
| Sterile | no |
| Re-usable | no |
| Patient contact | no |
| Contact with blood and blood plasma | yes |

Classification

Europe CE-IVD
CE labeling

In accordance with Directive 98/79/EC on in vitro diagnostic medical devices as accessory.

Article which, whilst not being an in vitro diagnostic medical device, is intended specifically by Ritter GmbH to be used together with a device to enable that device to be used in accordance with its intended purpose.

Purity




Produced under cleanroom conditions
(ISO 14644-1 class 8)

Materials

| | |
|------|--------------------|
| Tip | PP – Polypropylene |
| Tray | PC – Polycarbonate |


Tip Dimensions

| | |
|---------------|---------|
| Tip Length | 97.0 mm |
| Concentricity | ≤ 1.5 |



Tray Dimensions

| | |
|-------------|----------|
| Tray length | 126.9 mm |
| Tray width | 88.4 mm |



Packaging

| | |
|--------|-------------------------|
| Tray | 96 tips |
| Box | 100 trays (9,600 tips) |
| Pallet | 12 boxes (115,200 tips) |

Specifications are subject to change without prior notice.

In case of any questions, please contact our sales representative via phone, e-mail or our website.

<https://www.ritter-medical.de/en/contact>