



clearKnights, 50 µl, DP




clearKnights, 200 µl, DP

clearKnights

Robotikspitzen für Tecan und weitere OEM-Workstations
Robotic tips for Tecan and other OEM workstations

Artikel-Nr.	Artikelbezeichnung	Tray-Farbe	Typ	VE	Stück pro Umkarton
Article no.	Article description	Tray color	Type	PU	Pieces per cover. box
49106 - 0000	clearKnights, 50 µl, tray, DP	gelb yellow	Tecan	10 x 96	4800
49102 - 0000	clearKnights, 200 µl, tray, DP	blau blue	Tecan	10 x 96	4800
49103 - 0000	clearKnights, 1000 µl, tray, DP	rot red	Tecan	10 x 96	4800

Artikel-Nr.	Artikelbezeichnung	Tray-Farbe	Typ	VE	Stück pro Umkarton
Article no.	Article description	Tray color	Type	PU	Pieces per cover. box
49106 - 1102	clearKnights, 50 µl, pre-sterilized, with filter, tray, DPS	gelb yellow	Tecan	10 x 96	4800
49102 - 1102	clearKnights, 200 µl, pre-sterilized, with filter, tray, DPS	blau blue	Tecan	10 x 96	4800
49103 - 1002	clearKnights, 1000 µl, pre-sterilized, tray, DPS	rot red	Tecan	10 x 96	4800
49103 - 1102	clearKnights, 1000 µl, pre-sterilized, with filter, tray, DPS	rot red	Tecan	10 x 96	4800



50 µl 200 µl 1000 µl



Typ „Tecan“ ist kompatibel mit den folgenden Robotikinstrumenten:

Type "Tecan" is compatible with the following robotic instruments:

- Tecan: Genesis, Freedom EVO, MiniPrep, Cavro
- Qiagen/Corbett: CAS 1200/4200/4800, QIAgility
- Perkin Elmer (Packard): Multiprobe II, Janus (mit Multiprobe-Kopf | with Multiprobe head)