

**Tip 10µl HYBRIDRACK sterile**

**Part Number : 072120**




**Ultra micro and standard Tips**

Universal fit on research grade pipettors  
Unfiltered tips (for filter tips, please see MultiGuard™ Barrier Tips)  
Accurate, precise and reliable  
Exceptional clarity  
Extended length tips to reach bottom of tubes  
Certified as DNase, RNase, PCR inhibitor, DNA free

**HybridRack**

- Paperboard base easily recyclable
- 70% less plastic than a traditional tip rack → 8 lb / 3,6 kg of plastic saved per case of 300ul tips !
- Combination of sturdy polypropylene deck and strong paperboard base allows for stable pipetting directly from HybridRack, even with Multi-Channel
- Pack box is made from 100% recycled material, 50% post consumer content
- Plastic Deck is compatible with Sorenson Platinum and Sorenson MultiRack tip racks
- Color coded decks depending on tip volume

Volume : 10  $\mu$ l  
Shape : LongReach  
Color : Natural  
Filter : No  
Sterile : Yes  
Graduated : Graduated : 2, 5, 10  $\mu$ l  
Presentation : HybridRack  
Low binding : No  
Length : 45,34 mm  
Top diameter : 6,35 mm  
Bottom diameter : 0,38 mm



Sales unit : 10 racks of 96 tips

