


**Tip 300µl HYBRIDRACK 59mm sterile**

**Part Number : 072140**



**Ultra micro and standard Tips**

Universal fit on research grade pipettors  
Unfiltered tips (for filter tips, please see MultiGuard™ Barrier Tips)  
Accurate, precise and reliable  
Exceptional clarity  
Extended length tips to reach bottom of tubes  
Certified as DNase, RNase, PCR inhibitor, DNA free

**HybridRack**


- Paperboard base easily recyclable
- 70% less plastic than a traditional tip rack → 8 lb / 3,6 kg of plastic saved per case of 300ul tips !
- Combination of sturdy polypropylene deck and strong paperboard base allows for stable pipetting directly from HybridRack, even with Multi-Channel
- Pack box is made from 100% recycled material, 50% post consumer content
- Plastic Deck is compatible with Sorenson Platinum and Sorenson MultiRack tip racks
- Color coded decks depending on tip volume



**DOMINIQUE  
DUTSCHER  
DISTRIBUTION**

# Data sheet

Volume : 300  $\mu$ l  
Shape : NX  
Color : Natural  
Filter : No  
Sterile : Yes  
Graduated : Graduated : 10, 50, 100  $\mu$ l  
Presentation : HybridRack  
Low binding : No  
Length : 58,57 mm  
Top diameter : 7,49 mm  
Bottom diameter : 0,51 mm



Sales unit : 10 racks of 96 tips

