


Tip 300µl HYBRIDRACK XL 78 mm sterile

Part Number : 072150



Ultra micro and standard Tips

Universal fit on research grade pipettors
Unfiltered tips (for filter tips, please see MultiGuard™ Barrier Tips)
Accurate, precise and reliable
Exceptional clarity
Extended length tips to reach bottom of tubes
Certified as DNase, RNase, PCR inhibitor, DNA free

HybridRack

- Paperboard base easily recyclable
- 70% less plastic than a traditional tip rack → 8 lb / 3,6 kg of plastic saved per case of 300ul tips !
- Combination of sturdy polypropylene deck and strong paperboard base allows for stable pipetting directly from HybridRack, even with Multi-Channel
- Pack box is made from 100% recycled material, 50% post consumer content
- Plastic Deck is compatible with Sorenson Platinum and Sorenson MultiRack tip racks
- Color coded decks depending on tip volume



**DOMINIQUE
DUTSCHER
DISTRIBUTION**

Data sheet

Volume : 300 μ l
Shape : XT
Color : Natural
Filter : No
Sterile : Yes
Graduated : Graduated : 10, 50, 100, 200 μ l
Presentation : HybridRack
Low binding : No
Length : 77,50 mm
Top diameter : 7,6 mm
Bottom diameter : 0,46 mm

Sales unit : 8 racks of 96 tips

