

Tip 1000µl HYBRIDRACK sterile

Part Number : 072160




Ultra micro and standard Tips

Universal fit on research grade pipettors
Unfiltered tips (for filter tips, please see MultiGuard™ Barrier Tips)
Accurate, precise and reliable
Exceptional clarity
Extended length tips to reach bottom of tubes
Certified as DNase, RNase, PCR inhibitor, DNA free

HybridRack

- Paperboard base easily recyclable
- 70% less plastic than a traditional tip rack → 8 lb / 3,6 kg of plastic saved per case of 300ul tips !
- Combination of sturdy polypropylene deck and strong paperboard base allows for stable pipetting directly from HybridRack, even with Multi-Channel
- Pack box is made from 100% recycled material, 50% post consumer content
- Plastic Deck is compatible with Sorenson Platinum and Sorenson MultiRack tip racks
- Color coded decks depending on tip volume

Volume : 1250 μ l
Shape : XL
Color : Natural
Filter : No
Sterile : Yes
Graduated : Graduated : 200,500,1000
Presentation : HybridRack
Low binding : No
Length : 102,00 mm
Top diameter : 8,89 mm
Bottom diameter : 0,81 mm



Sales unit : 8 racks of 96 tips

