


PDS No. 77536x	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 01	Sapphire Pipette Tips, 200 µl, Universal, Low Retention				
	Greiner Item-No. 77536x				
Valid for Item-No.:	<b>775360</b>	<b>775361</b>	<b>775362 (sterile)</b>	<b>775364</b>	<b>775365 (sterile)</b>

1.	Description / Specification	
1.1	Description	PP Universal, Low Retention, Sapphire Pipette Tip, 200 µl, graduated For liquid handling of volumes in the range of 2 µl up to 200 µl, universal fit, graduated  775360: bag, non-sterile 775361: rack, non-sterile 775362: rack, sterile 775364: refill, non-sterile 775365: refill, sterile  Matching Rack: 770320 (yellow insert)
1.2	Dimensions	See Customer Drawing
1.3	Volume	2-200 µl
1.4	Material / Resin	Pipette tip: PP (low retention Polypropylene), free of heavy metal Rack/lid: PP (Polypropylene), free of heavy metal
1.5	Colour	Pipette tip: translucent Rack/lid: translucent; Rack insert: yellow
1.6	Sterilization	775360, -361, -364: No 775362, -365: by irradiation [SAL10 <sup>-6</sup> ]
1.7	Quality Control	- Raw Material control: physical testing - Product control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for pipetting of liquids to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Universal tip suitable for commonly used laboratory pipettors such as Eppendorf, Gilson, Sapphire, Brand - see valid compatibility list

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA, PCR inhibitors and endotoxins (pyrogens)
2.2	Temperature range	Room Temperature (15°C - 30°C)
2.3	Autoclavability	775360, -361, -364: 121°C / 20 min. 775362, -365: n.a. (sterile by irradiation)
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-sevices/download-center.html">https://www.gbo.com/en_INT/know-how-sevices/download-center.html</a>
2.6	Shelf life	3 years
2.7	Other Information	Compatibility testing procedure of the pipette tips follows the ISO 8655

3.	Packaging	775360	775361, 775362	775364, 775365
3.1	Pieces / packaging unit	1000 (bag)	96 (rack)	96 (rack / refill unit)
3.2	Pieces / Box	-	960	n.a.
3.3	Pieces / Case	15.000	5.760	4.800
3.3	Lot-No.	A YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)		
3.4	Other Information	Certificate of Quality to download		

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 00	Date 17 May 2021	Date 18 May 2021	Date 19 May 2021	
Date 20.11.2020	Name S. Kaelberer	Name J. Gaiser	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.