

Technical Data Sheet

1. Product Description

Product code: 0911R00034
Description: pH Indicator Paper pH 0-10

2. Overview

pH Indicator paper supplied in various presentations including books, strips, rolls and sheets. pH indicator papers are manufactured by impregnating high quality paper with indicator or mixed indicator solutions. The paper is dried, cut and packed as required.

3. Instructions for use

- a. Immerse the indicator paper in the test solution for 1 second.
- b. Remove and compare the color produced against the nearest matching color on the color chart for result

4. Indication range

pH 0 - 1 - 2 - 3 - 4 - 5 - 6 - 7 - 8 - 9 - 10

5. Shelf Life

3 years from date of manufacture

6. Product Information

H S Code: 38220000

7. Safety Information according to GHS

Storage class: 10 - 13 Other liquids and solids

8. Storage and Shipping information

Store at +15 °C to +25°C

9. Specifications

Application test: Conforms