

500 μ L, Model 750 SN SYR, Cemented NDL, Custom gauge, length, point style Specification Sheet (Part# 80808)



Product Specifications

Volume	500 μ L
Syringe Type	Microliter
Syringe Series	700
Termination	Luer Tip Cemented Needle Special (LTSN)
Needle	Cemented Needle
Custom Gauge	Custom Gauge
Custom Point Style	Custom Point Style
Custom Needle Length	0.75 inch (19.1 mm)
Std. Needle Dead Volume	N/A
Fluid Path	SST, borosilicate glass
Autoclavable	No
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.135 inches (3.26 mm)

Product Specifications

Barrel Outer Diameter	0.305 inches (7.75 mm)
Sales Quantity	1 syringe per pack

HAMILTON 

Web: www.hamiltoncompany.com

USA: +1-800-648-5950

Europe: +41 58 610 10 10

Hamilton Americas & Pacific Rim

Hamilton Company Inc.
4970 Energy Way
Reno, Nevada 89502 U.S.A.
Tel: +1-800-648-5950
Fax: +1-775-856-7259
sales@hamiltoncompany.com

Hamilton Europe, Asia, & Africa

Hamilton Bonaduz A.G.
Via Crusch 8
CH-7402 Bonaduz, Switzerland
Tel: +41 58 610 10 10
Fax: +41-58-610-00-10
contact@hamilton.ch

©2016 Hamilton Company. All Rights Reserved.

All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.

Revision: A / 2010-05-11 Printed in the U.S.A.

To find a representative in your area, please visit hamiltoncompany.com/contacts