

1 μ L, Model 7001 KH SYR, Knurled Hub NDL, 25 ga, 2.75 in, point style 2

Specification Sheet (Part# 80135)



Product Specifications

Volume	1 μ L
Syringe Type	Microliter
Syringe Series	7000
Termination	Knurled Hub (KH)
Needle	KH
Gauge	25 gauge
Point Style	2
Needle Length	2.75 inch (69.9 mm)
Needle Inner Diameter	0.0057 inches (0.145 mm)
Needle Outer Diameter	0.020 inches (0.515 mm)
Std. Needle Dead Volume	0 μ L
Fluid Path	SST, tungsten
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	6000 psig (413.7 bar)
Maximum Temperature	115°C / 239°F

Product Specifications

Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.0057 inches (0.145 mm)
Barrel Outer Diameter	0.308 inches (7.82 mm)
Sales Quantity	1 syringe per pack

MINIQUÉ DUTSCHER SAS

© 2016 Hamilton Company. All Rights Reserved.

All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.

Revision: A / 2010-05-11 Printed in the U.S.A.

HAMILTON 

Web: www.hamiltoncompany.com

USA: +1-800-648-5950

Europe: +41 58 610 10 10

Hamilton Americas & Pacific Rim

Hamilton Company Inc.
4970 Energy Way
Reno, Nevada 89502 U.S.A.
Tel: +1-800-648-5950
Fax: +1-775-856-7259
sales@hamiltoncompany.com

Hamilton Europe, Asia, & Africa

Hamilton Bonaduz A.G.
Via Crusch 8
CH-7402 Bonaduz, Switzerland
Tel: +41 58 610 10 10
Fax: +41-58-610-00-10
contact@hamilton.ch

To find a representative in your area, please visit hamiltoncompany.com/contacts