

# 5 $\mu$ L, Model 7105 KH SYR, Knurled Hub NDL, 24 ga, 2.75 in, point style 3

Specification Sheet (Part# 88000 )



## Product Specifications

Volume	5 $\mu$ L
Syringe Type	Microliter
Syringe Series	7000
Termination	Knurled Hub (KH)
Needle	Knurled Hub (KH)
Gauge	24 gauge
Point Style	3
Needle Length	2.75 inch (69.9 mm)
Needle Inner Diameter	0.0128 inches (0.325 mm)
Needle Outer Diameter	0.022 inches (0.566 mm)
Std. Needle Dead Volume	0 $\mu$ L
Fluid Path	SST, tungsten
Autoclavable	Yes
Gas Sterilizable	Yes
Pressure	6000 psig (413.7 bar)
Maximum Temperature	115°C / 239°F

## Product Specifications

Minimum Temperature	10°C / 50°F
Barrel Inner Diameter	0.0128 inches (0.325 mm)
Barrel Outer Diameter	0.308 inches (7.82 mm)
Sales Quantity	1 syringe per pack

**HAMILTON** 

Web: [www.hamiltoncompany.com](http://www.hamiltoncompany.com)

USA: +1-800-648-5950

Europe: +41 58 610 10 10

#### Hamilton Americas & Pacific Rim

Hamilton Company Inc.  
4970 Energy Way  
Reno, Nevada 89502 U.S.A.  
Tel: +1-800-648-5950  
Fax: +1-775-856-7259  
[sales@hamiltoncompany.com](mailto:sales@hamiltoncompany.com)

#### Hamilton Europe, Asia, & Africa

Hamilton Bonaduz A.G.  
Via Crusch 8  
CH-7402 Bonaduz, Switzerland  
Tel: +41 58 610 10 10  
Fax: +41-58-610-00-10  
[contact@hamilton.ch](mailto:contact@hamilton.ch)

© 2018 Hamilton Company. All Rights Reserved.

All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.

Revision: A / 2010-05-11 Printed in the U.S.A.

To find a representative in your area, please visit [hamiltoncompany.com/contacts](http://hamiltoncompany.com/contacts)