

# 30 gauge, Kel-F Hub Needle, custom length (0.375 to 12 in), point style 2, 3, or 4, 6/PK

Specification Sheet (Part/REF # 7750-21)

Kel-F hub (KF) needles can be used with Luer Tip (LT) and PTFE Luer Lock (TLL) syringes and connectors. Standard needle length is 2" / 51 mm. Kel-F hub (KF) needles are available in convenient six packs. Custom lengths and point styles are available. These needles are available from 10 - 31 gauge.

## Product Specifications

<b>Gauge</b>	30 gauge
<b>Point Style</b>	2, 3, or 4
<b>Needle Length</b>	Custom Needle Length
<b>Autoclavable</b>	Yes
<b>Gas Sterilizable</b>	Yes
<b>Sales Quantity</b>	6 needles per pack
<b>Hub</b>	Kel-F

Specifications are subject to change without notice

<https://www.hamiltoncompany.com/laboratory-products/needles/7750-21>

**Spec. Version B**