

# 23 gauge, Kel-F Hub NDL, 2 in, point style 2, 6/PK

Specification Sheet (Part #90123)



Kel-F hub (KF) needles can be used with Luer Tip (LT) and PTFE Luer Lock (TLL) syringes and connectors. Standard needle length is 2" / 51 mm. Kel-F hub (KF) needles are available in convenient six packs. Custom lengths and point styles are available.

## Product Specifications

<b>Gauge</b>	23 gauge
<b>Point Style</b>	2
<b>Autoclavable</b>	No
<b>Hub</b>	Kel-F
<b>Sales Quantity</b>	6 needles per pack
<b>Length</b>	2 inch (51 mm)