

TECHNICAL DATA SHEET

Part number	075003
Description	Complete Drummond Pipet-Aid XP Gravity pipettor



Description :

- ✓ For glass and plastic pipettes from 1 to 100 mL
- ✓ Reliability, robustness and ease of maintenance make it the most acclaimed instrument in its category
- ✓ Identical to the Pipet-Aid XP with an additional gravity function (distribution by gravity, without pressure)
- ✓ Power supply : Rechargeable battery 9 V
- ✓ Speed : 3 speeds (normal, slow, fast) + gravity
- ✓ Specificity : Filtering nose for culture

Pieces by sales :

1 Pipettor
