



# Drummond Scientific Company

500 PARKWAY, BROOMALL, PENNSYLVANIA 19008-4203

PHONE : 610 – 353 – 0200

FAX 610 – 353 – 6204

[www.drummondsci.com](http://www.drummondsci.com)

## Drummond Scientific 100ul Calibrated Pipets Specification Sheet

Volume 100ul (+/- 1%)

Catalog # 2-000-100

Color code: Blue

100ul Line Location: 2.953" (+/- .005")

Length: 5.000" (+/- .010")

OD: .0665" (+/- .0030")

ID: .05130" (+.00012/- .00015)

Glass type: Borosilicate

Packaging: 250 pieces per container

Aspirator Tube included

