



Drummond Scientific Company Technical Data Sheet

**Product # 1-000-0200/CA**  
**20.0ul Microcaps**

**Dimensions**

Volume: 20.00ul (+/- 1.00%)

Length: 2.520" +/- .005" (64mm (+/- .127mm))

Glass O.D. - .0355" +/- .0020" (.902mm (+/- .051mm))

Glass I.D. - .02483" +/- .00013" (.6307mm (+ .0030mm))

Glass type: Borosilicate

**Packaging**

Container: Printed Drummond Flex Tube

Vial Size: 3/4" dia. x 3 1/2" L

Foam Plug: 11/16" dia. x 15/16" L

Cap Size: 20mm w/ .056 hole

Quantity per vial: 100 glass pipettes.

Packaged with Bulb Assembly

10 vials per package

Material listed herein is produced in accordance with ASTM Standard Specification for Disposable Glass Micropipets number E 672-87. Standards used in manufacture are traceable to the National Institute of Standards and Technology.

Alexis Piccirilli  
Quality Supervisor