


Dialmatic Microdispenser

- Adjustable, direct reading plunger sets volume
- Delivery scale engraved on plunger shaft
- Fractional deliveries possible with premarked indexes
- Volumetric accuracy of $\pm 1\%$ (at full volume)
- Borosilicate glass bores with stainless steel plungers, (PTFE tips on all plungers 25 μL or larger)

Cat. No.	Total Capacity (μL)	Recommended Range (μL)	Delivery Increments (μL)
3-000-210	10	2.0 to 10	1
3-000-225	25	5 to 25	1
3-000-250	50	10 to 50	2
3-000-275	100	20 to 100	10



Replacement Parts Kits for Dialmatic Microdispensers

Cat. No.	Microdispenser Volume Range (μL)	Dialmatic Microdispenser Cat. No.	Kit Contents
3-000-010-K	2-10	3-000-210	2 plungers, 2 O-rings, 1 collet, 1 coil spring, 1 spring guide
3-000-025-K	5-25	3-000-225	1 plunger, 2 O-rings, 1 collet, 1 coil spring, 3 PTFE tips, 1 spring guide
3-000-050-K	10-50	3-000-250	1 plunger, 2 O-rings, 1 collet, 1 coil spring, 3 PTFE tips, 1 spring guide
3-000-075-K	20-100	3-000-275	1 plunger, 2 O-rings, 1 collet, 1 coil spring, 3 PTFE tips, 1 spring guide

Microdispenser Accessories and Replacement Parts

Bulk Replacement Parts

Cat. No.	Description	Pack
PTFE® Tips		
3-000-0101	PTFE Tips, 25 μL	50
3-000-0102	PTFE Tips, 50 μL	50
3-000-0103	PTFE Tips, 100 μL	50

NOTES: (1) PTFE tips are not required on 20 μL and smaller units.
 (2) Plunger wire and O-Ring are not required on 1000 μL size.

Replacement Bores for All Microdispensers

Size (μL)	Glass Color Code	Bores	Pack
1 to 5	White	3-000-105-G	100
10	Orange	3-000-210-G	100
20	Black	3-000-320-G	100
25	Double White	3-000-225-G	100
50	Green	3-000-250-G	100
100	Blue	3-000-275-G	100
100/200	Blue	3-000-275-G	100
250	Black	3-000-388-G	100
1000	Plastic Tip	3-000-590-G	100