



**Drummond Scientific Spec Sheet Catalog # 8-000-7520/5
Plasticrit Plastic 75 mm Hematocrit Tubes - Plain**

Length: 2.953“(± .020”)

Outside Diameter: .0740” (± .0030”)

Inside Diameter: .04300” (± .00200”)

Quantity per package: 5 vials/ box – 200 tubes per vial

Vial is labeled with Julian Lot number

Note: All components are manufactured according to the Drummond Quality Manual Standards using measuring equipment that is NIST Traceable.

DOMINIQUE DUTSCHER SAS