

10 μ L, Model 1701 N SYR Cemented NDL 26s ga 2 in point style 2

Specification Sheet (Part# 80000)



Product Specifications

Volume	10 μ L
Syringe Type	Gasight
Syringe Series	1700
Termination	Cemented Needle (N)
Gauge	26s gauge
Point Style	2
Needle Length	2 inch (51 mm)
Needle Inner Diameter	0.005 inches (0.13 mm)
Needle Outer Diameter	0.0187 inches (0.47 mm)
Std. Needle Dead Volume	0.64 μ L
Fluid Path	SST, borosilicate glass, PTFE
Autoclavable	No
Gas Sterilizable	Yes
Pressure	1000 psig (68.9 bar)
Maximum Temperature	50°C / 122°F
Minimum Temperature	10°C / 50°F

Product Specifications

Barrel Inner Diameter	0.018 inches (0.460 mm)
Barrel Outer Diameter	0.260 inches (6.60 mm)
Sales Quantity	1 syringe per pack

HAMILTON 

Web: www.hamiltoncompany.com

USA: +1-800-648-5950

Europe: +41 58 610 10 10

Hamilton Americas & Pacific Rim

Hamilton Company Inc.
4970 Energy Way
Reno, Nevada 89502 U.S.A.
Tel: +1-800-648-5950
Fax: +1-775-856-7259
sales@hamiltoncompany.com

Hamilton Europe, Asia, & Africa

Hamilton Bonaduz A.G.
Via Crusch 8
CH-7402 Bonaduz, Switzerland
Tel: +41 58 610 10 10
Fax: +41-58-610-00-10
contact@hamilton.ch

© 2018 Hamilton Company. All Rights Reserved.

All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.

Revision: A / 2010-05-11 Printed in the U.S.A.

To find a representative in your area, please visit hamiltoncompany.com/contacts