

FICHE TECHNIQUE

Référence	077503
Désignation	Pointe 10µL gel saver rack

Description :

Ultra Micro Gel Loading Tip

Volume: 0.5-10µl

For use with but not limited to:

Eppendorf®

Nichiryo

Oxford-Benchmate®

Socorex



Product No.	Color	Package Type	Qty/ Unit	Qty/ Case
124-R204	Natural	204 Tip Rack	1,020	10,200*
124-R204S	Natural	204 Tip Rack - S	1,020	10,200*

11

w - with
w/o - without

H - Hinged
S - Sterile

* - may be purchased by the unit

Unité de vente :

204 Pointes