

BEAD RUPTOR 96 WELL PLATE ADAPTER - SET OF TWO

Sku: 27-101



GENERAL INFO

Bead Ruptor 96 Well Plate Adapter - Set of Two

The Bead Ruptor 96 Well Plate Adapter Set contains two pairs of stainless steel plates. Each pair will hold 3 x standard well plates, 2 x 1 mL deep well plates or 1 x 2 mL deep well plates

A total of six standard well plates, 4 x 1 mL well plates, or 2 x 2 mL deep well plates can be processed in a single run using the well plate adapter set on the Bead Ruptor 96.

The Bead Ruptor 96 well plate adapter holds well plates or racks of strip well tubes securely during processing.
Ordering Information: Includes a pair of two stainless steel well plate adapters.

SPECIFICATIONS

Capacity:	3 x standard well plates, 2 x 1 mL deep well plates, 1 x 2 mL deep well plate
Volume Range:	200 µL-1.5mL
Compatibility:	Any SBS format well plate
Plate Sealing:	Recommended to seal plate with silicone map