



MICROTUBE OR SCREW CAP TUBE HOLDER

SKU: 27-106

GENERAL INFO

Bead Ruptor 96 PTFE Tube Adapter - Set of Two

The Bead Ruptor 96 Tube Adapter Set contains two pairs of PTFE tube adapters. Each pair will hold 24 x 0.5 mL screw cap tubes, 24 x 2 mL screw cap tubes, 24 x 1.5 mL microcentrifuge tubes or 24 x 2 mL microcentrifuge tubes.

A total of 48 tubes can be processed in a single run using the tube adapter set on the Bead Ruptor 96.

The Bead Ruptor 96 tube adapter holds microcentrifuge or screw cap tubes securely in place during processing.
Ordering Information: Includes a pair of two PTFE tube adapters.

SPECIFICATIONS

Capacity: 24 x 2 ml screw cap, 24 x 1.5 mL microcentrifuge tubes, 24 x 2 mL microcentrifuge tubes

Volume Range: 200 µL-1.5mL

Compatibility: Omni 0.5 mL, 1.5 mL or 2 mL pre-filled bead tubes