


2 ML DEEP WELL PLATES PRE-FILLED WITH 5 MM STAINLESS STEEL BEADS – 5 PACK

Sku: 27-6005



GENERAL INFO

Pre-filled 2 mL 96 deep well plates with 5 mm metal beads

Bead Ruptor bead beating tubes are pre-filled with an optimized amount of lysing matrix to provide efficient and high yield sample homogenates.

Product Features:

- 2.0 mL polypropylene 96 deep well plates. Wells pre-filled with 5.0 mm metal beads
- Recommended for: Operation on Bead Ruptor 96
- Certification: RNase & DNase Free

Ordering Information: Unit includes 5 pre-filled well plates sealed with silicon mats. 5 additional sealing mats included.

SPECIFICATIONS

Tube material:	Polypropylene.
Recommended for:	Processing samples on Bead Ruptor 96
Well volume:	2.0 mL
Bead material:	5.0 mm metal
Processing range :	0.2 mL to 1.5 mL
Temperature Range:	-80 to 121°C
Force:	4,000 x g
Certification:	RNase & DNase free