

## UNIVERSAL MICROBIAL GRINDING MIX (2 ML TUBES) NUCLEASE FREE - 50 PACK

SKU: 19-631



### GENERAL INFO

#### Universal Microbial Grinding Mix

Bead Ruptor bead beating tubes are pre-filled with an optimized amount of lysing matrix to provide efficient and high yield sample homogenates.

#### Product Features:

- 2.0 mL Polypropylene screw cap
- Tubes pre-filled with 0.5 mm ceramic beads
- Caps contain silicone O-ring

**Recommended for:** Larger micro-organisms: Protazoans, Nematodes, Parasites

**Certification:** RNase & DNase Free

**Ordering Information:** Unit includes 50 bead beating tubes prefilled with 0.5 mm ceramic beads

### SPECIFICATIONS

**Tube material:** Polypropylene

**O-ring material:** Silicone

**Tube volume:** 2.0 mL

**Bead material:** 0.5 mm ceramic

**Processing range :** 0.2 mL to 1.5 mL

**Temperature Range:** -80 to 121°C

**Force:** 20,000 x g

**Certification:** RNase & DNase free