

UNIVERSAL MICROBIAL GRINDING MIX (2 ML TUBES) NUCLEASE & MICROBIAL DNA FREE - 50 PACK

Sku: 19-631D



GENERAL INFO

Universal Microbial Grinding Mix

Microbial DNA free Bead Ruptor bead beating tubes are molecular biology grade, nuclease free tubes, pre-filled with an optimized amount of sterilized lysing matrix to provide efficient and high yield sample homogenates. Products are rigorously tested using microbial DNA PCR assays.

Product Features:

- 2.0 mL Polypropylene screw cap tubes pre-filled with 0.5 mm ceramic beads
- Ideal for applications that are sensitive to the presence of bacterial or fungal DNA contamination.
- Recommended for: Larger micro-organisms: Protozoans, Nematodes, Parasites
- Certification: RNase & DNase Free, Microbial DNA Free

Ordering Information: Unit includes 50 bead beating tubes prefilled with 0.5 mm ceramic beads

- **Certification:** RNase & DNase free, Microbial DNA free

SPECIFICATIONS

Tube material:	Polypropylene
O-ring material:	Silicone
Tube volume:	2.0 mL
Bead material:	0.5 mm ceramic
Processing range :	0.2 mL to 1.5 mL
Temperature Range:	-80 to 121°C
Force:	20,000 x g
Certification:	RNase & DNase free, Microbial DNA Free