




## REPLACEMENT PARTS

1-888-850-6208

[www.scientificindustries.com](http://www.scientificindustries.com)

### **0K-0400-901 KIT FOR PCR PLATES FOR MICROPLATE GENIE AND MULTI- MICROPLATE GENIE**

Two black, cylindrical bumpers are shown against a white background. They are positioned vertically, one slightly to the left and one slightly to the right of the center. They appear to be made of a solid material, possibly plastic or metal, and have a smooth, slightly reflective surface. The lighting creates soft shadows to the left of each bumper.

These longer, stiffer bumpers replace the existing bumpers on the MicroPlate Genie, Digital MicroPlate Genie Pulse Multi-MicroPlate Genie, and Digital Multi-MicroPlate Genie Pulse Mixers, to accommodate the taller PCR Plates. Please note that (4) kits are necessary for Multi-MicroPlate Genie and Digital Multi-MicroPlate Genie Pulse. If the serial number on your MicroPlate Genie is lower than M40-2075 you will also need to purchase Part No. 0K-0400-900 Shock Absorbing Feet 4 pk. If the serial number on your Multi-MicroPlate Genie is lower than M4K-1390 you will also need to purchase Part No. 0K-4000-900 Shock Absorbing Feet 4 pk.

Print this Page