

**Product Name:** **Vacucap™ Tube Closures**

**Catalogue No.:** **T402-13 Series**

This document replaces any previous version

## 1. Product Description:

- **Vacucap™ Tube Closure for 13 mm Tube:** Non sterile, and disposable.
- **Colors:**
  - Gray



## 2. Packaging:

- **Case:** 1 000 units in a package / 6 000 units per case

## 3. Product Specifications:

- Material: Polyethylene;
- Temperature range: -80°C to +80°C;
- Not autoclavable;
- Fit on 13 mm glass test tubes and plastic culture tubes;

## 4. Standards and Conformity:

- **ISO 2859-1:** Sampling and inspection procedures
- **FDA:** Resin conforms to FDA 21 CFR 177.1520
- **CONEG:** Plastics and colorants are in conformity with CONEG standards for heavy metals

## 5. Quality Assurance:

- Clear, no presence of contamination in plastic;
- Visual attributes;
- Dimension measurements;
- Closure verification.

## 6. Traceability:

- **Lot No. Composition:** 8 to 9 digits
- **The lot number can be found in one or all of these locations:**
  1. On exterior case label;
  2. On label inserted inside the master case;
  3. On the inner bag.

## 7. Storage Conditions:

- Store at room temperature in normal warehouse conditions;
- Avoid temperature variations and humidity;
- Protect from any possible contamination;
- Protect from any damage to the packaging.

## 8. Recommended Use:

- Follow chemical resistance chart recommendations.

### **OWNERSHIP OF MATERIALS:**

*Materials and information contained in this document are **Simport Scientific's** copyrighted and are protected by worldwide copyright laws and treaty provisions. They may not be copied, reproduced, modified, published, uploaded, posted, transmitted or distributed in any way.*