

# TECHNICAL DATA SHEET

## EUROTUBO® 12 ml screw cap PE tubes, round bottom

codes 301402 - 301403

code	description	sterile	maximum volume (ml)	rec. volume (ml)	case quantity	case weight (kg)	case volume (m³)	cases pallet
301402	polystyrene individually wrapped	STERILE EO	14.4	12.0	6 x 250	13.40	0.096	16
301403	polystyrene	NO	14.4	12.0	6 x 250	12.90	0.082	20

Medical device for in vitro diagnostic, Directive 98/79/EC




### SPECIFIC USE

Made of polystyrene. Green cap made of high density polyethylene.  
 Dimensions: 15 x 102 mm. The external skirt allows the tubes to remain free-standing. Supplied screwed. Recommended volume: 12 ml. The sterile model (ethylene oxide) is supplied individually in flow-pack bag, with indication of batch, expiration date, etc.  
 According to the guidelines for sterile products. Resistance to centrifugation: 7,500 xg.

### MEASUREMENTS

external cap Ø mm A	external tube Ø mm B	internal tube Ø mm C	length with cap mm D	length without cap mm E
20.9	16.3	14.4	103.9	102.5

Measurements in mm (±0.09)



WRITTEN



David Campos  
Marketing Department

VERIFIED



Anna Mir  
Technical, Development and Quality Director

AUTHORISED



Mònica Torras  
Sales & Marketing Director