

## Technical Data Sheet

**Product:** Steril stopper C 40  
**Code-no.:** 890 57 40

---

**Material:** Silicosen®

**Quantity of airflow:** 2.000 - 4.500 ml/ min

**Dimension:** Length: 28 mm  
Diameter d1: 36 mm  
Diameter d2: 27 mm  
Diameter d3: 40 mm

**NS:** 35/ 45

**Evaporation of Moisture:** 0,6 ml in 10 days at 30°C

**Sterilization:** Dry air sterilization at 180°C for 20 - 30 minutes  
and high-pressure steam at 121 °C for  
20 minutes are possible.

The number of sterilizations is not limited.

Hirschmann Laborgeräte  
GmbH & Co. KG



.....  
Roland Benz  
Quality Manager

Eberstadt, 21.02.2020