

Technical Data Sheet

Product: Steril stopper S 28
Code-no.: 890 58 28

Material: Silicosen®

Quantity of airflow: 30 - 120 ml/ min

Dimension: Length l1: 18 mm
Length l2: 27 mm
Diameter d1: 28 mm
Diameter d2: 17 mm

NS: 17 - 26

Evaporation of Moisture: 0,55 ml in 10 days at 30°C

Sterilisation: Dry sterilized (180°C, for 15 - 30 minutes) and high-pressure steam sterilized

Washing: The stopper can be washed with hot water or with neutral detergent.

Hirschmann Laborgeräte
GmbH & Co. KG



.....
Roland Benz
Quality Manager

Eberstadt, 25.07.2018