










ACC-FL180180NL-11

	Reference Cryopal	ACC-FL180180NL-11
	Reference SAP	116523
	Designation	Flexible hose DN10 - 180/180 for Nitrogen Lg=3m equipped at both ends with a 180° wing nut.
	Dimensions / Size	1,1 m
	Material	Stainless steel Flexible metal "HYDRA" piping standard version, corrugated with stainless steel parallel corrugations reinforced by stainless steel braid.
	Use	Transfer of liquid nitrogen Max. operating pressure 32 bar Maximum operating temperature: -196 °C
	Origin	France
	Warranty	2 years
	Packing	Plastic sheath
	Weight	0,950 kg
	Regulation - Standard	Hydraulic test at 48 bar Conformity NF L 00-015C
	Benefits	Both ends: Bronze wing nut 3/4" - 16 UNF



DOMINIQUE DUTSCHER SAS