

## PRODUCT DATA SHEET

### Cryogenic gloves

Cryopal Commercial Reference	ACC-SECU-15 / ACC-SECU-16 / ACC-SECU-17 / ACC-SECU-18
SAP internal reference	151015 / 151016 / 151017 / 151018
Designation / Description	Cryo gloves - (Category III personal protective equipment)
Material	Glove and reinforcement thumb / index in water-repellent blue wool grain 1mm, lined in wool fleece, Hydrophilic and microporous insert. Seams made of polyamide thread. Cuff of 15cm in crust of cattle.
Origin	France
Packing	Plastic sheath
Use	Protective glove exclusively against mechanical hazards, industrial refrigeration and liquefied gas projections. No prolonged contact with liquefied gases.
Size	8 - 9 - 10 - 11
Regulation - Standard	EEC / 89/686 - EC certificate N ° 0072/011/162/05/97/0139 EN 420: General requirements EN 388: Mechanical risks Abrasion: 2/4 Slicing cut: 1/5 Tear: 2/4 Perforation: 2/4 EN 511: cold Convective cold: 2/4 Cold contact: 2/4 Impermeability to water: 1/1

