








ACC-RACK-5_RACK VERTICAL 12 LEVELS FOR TUBES 1,2/2ml

Version 0

	Benefits	<p>Made of aluminium, the rack combines lightness and good thermal conduction to promote the diffusion of refrigeration.</p> <p>The individual closures for each level allow the removal of a box of cryo-tubes without releasing the other boxes.</p>
	Cryopal Reference	ACC-RACK-5
	SAP Reference	116691
	Designation	Rack vertical 12 levels 1.2/2ml tube (for 81 or 100 boxes)
	Dimensions / Size	<p>Total height = 674 mm Width = 144 mm Depth = 136 mm Level height = 56 mm</p>
	Material	Aluminium
	Use	Store 1.2/2ml tubes in 81 or 100 cryo-boxes.
	Capacité de stockage	<p>12 cryo boxes of 81 or 100 tubes 1,2/2ml With 81 boxes = 972 tubes With 100 boxes = 1200 tubes</p>
	Warranty	1 year

	Origin	France
	Packaging	Plastic packaging
	Weight	1,650 kg
	HS Code	76169990
	Packing	0,675 x 0,144 x 0,137 m
	Regulation - Standard	Medical device class 1
	Maintenance Cleaning	-

For any information, feel free to contact us :



Cryopal
 Parc Gustave Eiffel - 8, avenue Gutenberg
 CS10 172 - Bussy Saint-Georges
 77 607 - Marne la Vallée cedex 3 - France



00 33 1 64 76 15 00



contact.cryopal@cryopal.com



<http://www.cryopal.com>