



**SUPER STRENGTH RED  
BIOHAZARD DISPOSAL BAGS WITH  
WARNING LABEL/STERILIZATION  
INDICATOR; 2.0MIL THICK, 2-4  
GALLON CAPACITY, POLYPROPYLENE  
(PACK OF 200)**

Extra Thick, Extra Strong High Temperature Bags  
Reduce Autoclave Time



Catalog No. : F13165-1419  
UPC No. : 191634001782

Protect against dangerous tears and leaks with extra strong, high impact polypropylene biohazard disposal bags.

- Red bags printed with standard biohazard warning symbol and precautionary procedures in four languages: English, French, German and Spanish
- Pass ASTM 1709 Dart Impact and ASTM 1922 Tear Resistance Tests
- Extra heavy-duty; 0.05mm (2 mil) thick blended polymer
- Autoclavable to 135°C (275°F)
- Sterilization indicator patch darkens when exposed to steam sterilization for 20 minutes
- Extra large sizes available, call for quote

Catalog No.:	Thickness	Dimensions (W X H)	Gallon Capacity	Quantity / Pack
F13165-1419	0.05mm (2.0mil)	36 x 48cm (14 x 19")	2 to 4	200