



Wayne, New Jersey 07470
Tel: 973-694-0500 * Fax: 973-694-7199
Internet Address: <http://www.bel-art.com>

Technical Data Sheet

December 11, 2012

Product #: 36940-0000

Description: "Sterileware" Sterile Sampler Spoons

Material: High Impact Polystyrene, sterilized in individual polyethylene bags. White polystyrene Colorant, materials meet FDA regulations 21 CFR 177.1640, 178.3297, 175.300

Approximate Dimensions: Overall 7" (178mm), Handle 6.12", Bowl Dimensions Width 0.81" (20mm), Length 0.96" (24.4mm), 0.26" (6.60mm) Deep, Volume approx-1.23ml

Other: No testing is done for Pyrogens / End toxins at this time
Based on the process and materials involved in the manufacturing of this product, this product is considered to be BSE/TSE protein free,
Based on the information provided to us this product does not contain or come in contact with lead, latex, Mercury, Melamine, formaldehyde or other contaminants
To the best of our knowledge based on information provided to us this product does not contain or come in contact with any animal origin.
Product should remain in its irradiated state as long as outer packaging remains completely intact. ANSI/AAMMI/ISO 11137-1994 Method 1, guidelines are followed.

Manufactured: Bel-Art Products , 2024 Broad Street, Pocomoke City MD, USA

*Carl Sisco
Bel-Art Products
Technical Support Mgr.
Quality Assurance
800-423-5278 x 4269
Fax 800-545-9796
csisco@belart.com*