



Bel-Art Products, a part of SP
 661 Route 23 South, Wayne, New Jersey 07470
 Tel: 973-694-0500 * Fax: 973-694-7199
www.sp-belart.com

DATA SHEET

This is to verify that all items noted below are in conformance with the Purchase Order Contract, drawings, specifications, and other applicable documentation; that all process certifications, material certifications, are on file at this or another Bel-Art facility and are available for review by the customer and or Government Representatives upon request.

Product #: 36930-0000 **Mfg. Date:** n/a

Description: "Sterileware" Sterile Sampler Spatula


Material: High-Impact Polystyrene Natural, White Polystyrene Colorant sterilized in individual polyethylene bags Scoop Material meets FDA requirements for regulation 21 part 177.1640. White polystyrene Colorant meets FDA regulations 21 CFR 177.1640, 178.3297, 175.300

Dimensions: Overall 229mm Spatula only Length-152mm, Width-19mm, Depth 13mm Volume approx-15cc.

Other: Irradiation exposure dose (avg.)= 23.5- 50.0 kGy. No testing is done for Pyrogens/ End toxins at this time. Based on the process and materials involved in the manufacturing of this product, this product is considered to be BSE/TSE protein free. AAMI guidelines are followed. Product should remain in its irradiated state as long as outer packaging remains completely intact. Product should be stored in room temperature, dry location. Based on the information provided to us this product does not contain or come in contact with lead, latex, Mercury, Melamine, formaldehyde any or other contaminants. To the best of our knowledge based on the information provided to us this product does not contain or come in contact with any animal origin.

Manufactured by: Bel-Art Products, a part of SP
2024 Broad Street Pocomoke City, MD 21851
An ISO-9001:2015 facility (certificate US19/819943409)
Made in the USA

Reviewed and Approved By:

Signature: 

Printed Name: Michael Grzelak

Title: Technical Support Analyst

Date: 3/26/2021