

## CERTIFICATE OF ANALYSIS

This is to verify that all items noted below are in conformance with the Purchase Order Contract, drawings, specifications, and other applicable documentation; that all process certifications, material certifications, are on file at this or another Bel-Art facility and are available for review by the customer and or Government Representatives upon request.

**Product #:** 36725-0000      **Mfg. Date:** NA      **Lot #:** General

**Description:** Sampler Spoon

**Material:** Polypropylene, Natural, 20% Talc Filled

**Dimensions:** Overall 7" (178mm), Handle 5.59" (139.75mm) Bowl 1.18" (29.5mm) wide x 1.47" (37.3mm) long x 0.36 (9mm) deep, 4.93ml Capacity

**Other:** Steam Autoclavable @ 121°C for 20 minutes. Materials meet FDA regulation 21 CFR 177.1570, Based on the process and materials involved in the manufacturing of this product, this product is considered to be BSE/TSE protein free. To the best of our knowledge based on information supplied to us by the resin manufacturers this product does not contain or come in contact with any animal origin.

Manufactured: Bel-Art Products Pocomoke City MD, USA

\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

Reviewed and Approved By:

**Signature:** *Carl Sisco*

**Printed Name:** Carl Sisco

**Title:** Quality Specialist

**Date:** 11/12/2015