



BEL-ART PRODUCTS

Pequannock, New Jersey 07440-1992 USA

Tel: 973-694-0500 Fax: 973-694-7199

Internet Address: <http://www.bel-art.com>

Data Sheet

Product #: 18905-0001

Description: Microcentrifuge Tube Ice Rack/Tray / Spill Containment Tray

Material: Polypropylene

Approximate Dimensions: Overall Outside Dimensions 9 1/8" L x 5 1/4"W x 2 7/8" H

Inside: Top 7.25" x 4.25" x 2.3125"D

Inside Bottom: 6.75" x 3.75"

2.08" Deep from top of grid to bottom of tray

Grid: ID 7.02" x 4.125" x .64" Deep from top of rack to top of grid

Hole Diameter in grid .47" (12mm)

Other: Grid consists of 5 rows x 10 rows = 50 places

Cover available (18905-0011)

Holds tube with less than 12mm Diameter / 1.5ml tubes