



BEL-ART PRODUCTS

Wayne, New Jersey 07470

Tel: 973-694-0500 * Fax: 973-694-7199

Internet Address: <http://www.bel-art.com>

Technical Data Sheet

Product #: 18796-0000

Description: 15 & 50ml Tube Holder

Material: Polypropylene 20% TALC Filled White

Approximate Dimensions: 3" OD x 2.65" H
15ml ID .64" (16.25mm)
50ml ID 1.24" (31.5mm)
.49oz / 13.7 grams

Other: Steam Autoclavable @ 121°C for 20 min. @ 15 psi

DOMINIQUE DUTSCHER SAS