



**Filter tip 20µl HYBRIDRACK low retention sterile**

**Part Number : 087250T**



**Multiguard barrier tips**

Sterilized Filtered tips, natural color

Ideal for PCR, DNA Amplification and radioisotope handling

Inert, the hydrophobic barrier prevents pipettor contamination

Universal fit on research-grade pipettors

Certified as DNase, RNase, PCR inhibitor and ADN free

**HybridRack**

- Paperboard base easily recyclable
- 70% less plastic than a traditional tip rack → 8 lb / 3,6 kg of plastic saved per case of 300ul tips !
- Combination of sturdy polypropylene deck and strong paperboard base allows for stable pipetting directly from HybridRack, even with Multi-Channel
- Pack box is made from 100% recycled material, 50% post consumer content
- Plastic Deck is compatible with Sorenson Platinum and Sorenson MultiRack tip racks
- Color coded decks depending on tip volume



**DOMINIQUE  
DUTSCHER  
DISTRIBUTION**

# Data sheet

Volume : 20  $\mu$ l  
Shape : Z tip  
Color : Natural  
Filter : Yes  
Sterile : Yes  
Gratuated : 10, 50, 100  $\mu$ l  
Presentation : HybridRack  
Low binding : Yes  
Length : 51,44 mm  
Top diameter : 7,42 mm  
Bottom diameter : 0,46 mm

Sales unit : 10 racks of 96 tips

20  $\mu$ l Standardization

