



Filter tip 300µl HYBRIDRACK low retention XL 78mm sterile

Part Number : 087280T



Multiguard barrier tips

Sterilized Filtered tips, natural color

Ideal for PCR, DNA Amplification and radioisotope handling

Inert, the hydrophobic barrier prevents pipettor contamination

Universal fit on research-grade pipettors

Certified as DNase, RNase, PCR inhibitor and ADN free

HybridRack


- Paperboard base easily recyclable
- 70% less plastic than a traditional tip rack → 8 lb / 3,6 kg of plastic saved per case of 300ul tips !
- Combination of sturdy polypropylene deck and strong paperboard base allows for stable pipetting directly from HybridRack, even with Multi-Channel
- Pack box is made from 100% recycled material, 50% post consumer content
- Plastic Deck is compatible with Sorenson Platinum and Sorenson MultiRack tip racks
- Color coded decks depending on tip volume



**DOMINIQUE
DUTSCHER
DISTRIBUTION**

Data sheet

Volume : 300 μ l
Shape : XT
Color : Natural
Filter : Yes
Sterile : Yes
Gratuated : 10, 50, 100, 200 μ l
Presentation : HybridRack
Low binding : Yes
Length : 77,5 mm
Top diameter : 7,6 mm
Bottom diameter : 0,46 mm



Sales unit : 8 racks of 96 tips

