


Filter tip 1250µl HYBRIDRACK low retention sterile

Part Number : 087290T



Multiguard barrier tips

Sterilized Filtered tips, natural color

Ideal for PCR, DNA Amplification and radioisotope handling Inert, the hydrophobic barrier prevents pipettor contamination Universal fit on research-grade pipettors

Certified as DNase, RNase, PCR inhibitor, ADN free and non pyrogenic.

HybridRack


- Paperboard base easily recyclable
- 70% less plastic than a traditional tip rack → 8 lb / 3,6 kg of plastic saved per case of 300ul tips !
- Combination of sturdy polypropylene deck and strong paperboard base allows for stable pipetting directly from HybridRack, even with Multi-Channel
- Pack box is made from 100% recycled material, 50% post consumer content
- Plastic Deck is compatible with Sorenson Platinum and Sorenson MultiRack tip racks
- Color coded decks depending on tip volume



**DOMINIQUE
DUTSCHER
DISTRIBUTION**

Data sheet

Volume : 1250 μ l
Shape : XL
Color : Natural
Filter : Yes
Sterile : Yes
Gratuated : 200,500,1000
Presentation : HybridRack
Low binding : Yes
Lenght : 102,00 mm
Top diameter : 8,89 mm
Bottom diameter : 0,81 mm



Sales unit : 8 racks of 96 tips

